


# 2MHz

## NMR Rock Core Analyzer




[www.magritek.com/rca](http://www.magritek.com/rca)

# Magritek NMR Rock Core Analyzer High Resolution and Fast

**NEW!**

## New probe for shales and tight rocks with 60 $\mu$ s TE

The new P34 high resolution 60  $\mu$ s TE (echo time) probe provides measurement capability not previously available on 2 MHz NMR rock core analyzers. The P34 probe is ideal for measurements of shales and tight rocks where critical information lies in  $T_2$  values below 1 ms.




## Measurements in minutes

Magritek's industry leading Q-switch probes have a very high signal-to-noise ratios, enabling measurements to be made in a very short time. Even low porosity samples can now be analyzed in a matter of minutes.

Mancos Shale					
Probe	TE	Acquisition time	# scans	SNR	Fluid Vol.
P49	100 $\mu$ s	1 minute	56	206	3.8 ml


## Measure cores with a wide range of diameters in the same analyzer

The Magritek NMR Rock Core Analyzer accepts a wide range of interchangeable probes, fitting sample diameters from under 1" up to full 4" cores. Additional options include overburden systems for high pressure and temperature, and strong, fast gradients for advanced measurements.



# Measure cores with a wide range of diameters in the same analyzer

The Magritek NMR Rock Core Analyzer incorporates a modern Halbach magnet providing a large clear bore (~120 mm) with excellent field uniformity. This large clear bore means that cores with a wide range of diameters can be measured in the same analyzer by using the five interchangeable probes listed below. With the additional option of overburden or gradients the analyzer is highly flexible. The fully insulated and temperature controlled magnet makes the analyzer very stable over a range of sample or environmental conditions.




## P34 probe: high resolution

Our new high resolution probe with 60  $\mu\text{s}$  TE (30  $\mu\text{s}$  tau) is ideal for shales and tight rock cores. The 60  $\mu\text{s}$  TE enables improved characterization for  $T_2$  times below 1 ms. It fits up to 1" diameter samples.

## P49, P54 and P71 probes: fast measurements

Our Q-switched series of probes provide superior signal-to-noise ratios enabling fast measurements. Magritek introduced Q-switch probe technology in 2009 setting a new standard for fast measurement probes. These probes enable measurements to be taken in minutes even with low porosity samples.



## P108 probe: full 4" core

Remarkable but true - with the new P108 probe the same analyzer can now measure full 4" (102 mm) diameter cores up to 100 mm in length.



Application	Probe	Maximum Sample diameter	Gradient	Overburden
Shales and tight cores	P34	1"	Yes	No
Standard cores	P49	1.5"		
Advanced measurements with overburden & gradients	P54	50 mm		
	P71	60 mm	No	Yes
Full 4" cores	P108	4"		

## Optional gradient

Magritek's gradients are strong (500 mT/m), fast switching, and have very low eddy currents providing the best platform for advanced diffusion and two dimensional measurement methods.

## Optional overburden and temperature control

The Magritek NMR Rock Core Analyzer works with most third party overburden systems using one of the standard probes. Custom probes are also available for specialized systems.



# Fully integrated software with every system




Magritek's hardware and software are seamlessly integrated - it just works! Magritek's **Dual Mode software** provides the optimal interface for your task. For standard measurements, we provide an interface making routine measurements straightforward. For the research user, the advanced mode provides complete flexibility with control of all parameters including the ability to write your own pulse sequences. The operation is transparent and efficient in both modes.

## Standard mode


- Stable parameter settings
- Simple interface
- Measurement queuing
- Easy to use

Project information and parameters



Run measurements and analyze data

All data easily available in Excel compatible format



## Advanced mode

For research applications, the advanced interface provides full control of all experimental parameters and enables the user to write their own pulse sequences for the Magritek NMR Rock Core Analyzer. Make use of an extensive library of built in functions or build your own using simple macro scripting. You can also fully customize the display and interface.

- Write your own pulse sequences
- Write your own macros
- Full display and interface control
- Powerful analysis features
- Full experimental control of all parameters

```

Editor (S) CPMGInt_pp_macropulst_program0
File Edit View Procedures Help

# Pulse sequence
intpp(d1r) # Reset internal parameter list



gradon(x,m3) # Set x shim
gradon(y,m4) # Set y shim
gradon(z,m5) # Set z shim

cleardata(m1)
pulse(mode,a1,p1,d1) # 90 RF pulse
delay(d2) # 90 - 180 delay
loop(l1,m1)
  pulse(mode,a2,p2,d1) # 180 RF pulse
  delay(d3) # 180 - acq delay
  acquire("integrate",m2,d4) # Acquire echo and wait
endloop("l1")

lst = endpp() # Return parameter list

# Phase cycle
phaseList = [2,0,2,0] # -x,x,-x,x : 90 phase
            [3,1,1,3] # -y,y,-y,y : 180 phase
            [2,0,0,0] # -x,x,-x,x : Acquire phase
    
```

Pulse programming is based around a simple scripting language which allows for the rapid development of pulse sequences as well as custom user interfaces.





## Technical Support & Service

### Service and Support

Magritek is committed to providing a range of superior service and support options to meet the different needs of our customers. Contact Magritek for details.

### Specifications

Operating frequency (Proton)	2 MHz
Maximum sample length	100 mm
Maximum sample diameter	102 mm (4")

### Standard system components

NMR spectrometer  
NMR rock core magnet  
Temperature controller  
Rock Core Analyzer software

### Standard weights and dimensions

Magnet weight	50 kg (110 lbs)
Magnet height	525 mm (20.6")
Magnet diameter	331 mm (13")
Spectrometer weight	5 kg (11 lbs)
Spectrometer dimensions	360 x 260 x 160 mm (14" x 10" x 2")
Temperature controller weight	3 kg (6.6 lbs)
Temperature controller dimensions	150 x 190 x 270 mm (6" x 7.4" x 10.6")

For further information please contact:

[www.magritek.com](http://www.magritek.com)

Magritek Limited  
32 Salamanca Road  
Wellington 6012  
New Zealand  
[sales@magritek.com](mailto:sales@magritek.com)

US: +1 (415) 287 0727  
UK: +44 (0) 20 328 7670  
HK: +852 8176 7196  
NZ: +64 4 920 7671  
Fax: +64 4 471 4665



**magritek**  
[www.magritek.com](http://www.magritek.com)