


## Measuring Lipophilicity with NMR



---

## Contents

Objectives	1
Introduction	1
Determining log $P_{ow}$ for some common solvents	2
Example $^1H$ spectra	3
References	9

## Objectives

The aim of this experiment is to determine the 1-octanol/water partition coefficient ( $P_{OW}$ ) of some common solvents (which will be the analyte in this experiment). Students will quantify the amount of the chosen analyte in water, before and after the addition of 1-octanol layer by NMR spectroscopy using the benchtop Spinsolve NMR spectrometer (Figure 1). This information will allow the calculation of the  $P_{OW}$  and  $\log P_{OW}$  of the analytes.

## Introduction

The lipophilic (“fat-loving”), or conversely the hydrophilic (“water-loving”) behaviour of a molecule is an important physicochemical property. From a pharmaceutical perspective, this property for a drug-like molecule strongly influences its distribution within the body. The typical method for determining the lipophilicity of a substance is to measure its distribution between two immiscible solvents, usually 1-octanol and water.

This allows for the calculation of the 1-octanol/water partition coefficient ( $P_{OW}$ ) and its logarithm ( $\log P_{OW}$ ). As stated by Lipinski’s “Rule of 5” for orally active drugs, a substance will have better drug-likeness if its  $\log P_{OW}$  is no greater than 5.<sup>1</sup> In this experiment, the  $\log P_{OW}$  for some common solvents will be experimentally determined by  $^1\text{H}$  NMR and compared to literature values.




Figure 1. Representative diagram of analyte distribution between the water and 1-octanol layers upon addition of 1-octanol.

## Determining log $P_{ow}$ for some common solvents

### Safety

Acetone, acetonitrile, 2-propanol, THF, ethanol, *tert*-butanol and 1-octanol are highly flammable, handle with caution.

### Procedure

Place approximately 200  $\mu\text{L}$  of acetone and **exactly** 500  $\mu\text{L}$  of water into a clean, dry NMR tube, cap the tube and invert to mix.

Record  $^1\text{H}$  NMR spectrum of the solution, obtain integration values for the acetone and water peaks. Normalise the water peak to a value 1.00.

Add **exactly** 500  $\mu\text{L}$  of 1-octanol to the NMR tube, cap the tube and invert 15-20 times to mix well. Allow the layers to separate (approximately 5 minutes).

Record another  $^1\text{H}$  NMR spectrum of the solution, and again obtain integration values for the acetone and water peaks. Normalise the water peak to a value 1.00.

Use integration values to calculate the  $P_{ow}$  and the log  $P_{ow}$  of acetone and write on a large whiteboard.

Using at least 4 other values (obtained by other students), calculate the mean and standard deviation.

Repeat the above procedure, replacing acetone with **one** of the following: ethanol, 2-propanol, acetonitrile, THF or *tert*-butanol.

### Calculations

The octanol/water partition coefficient ( $P_{ow}$ ) is equal to the concentration of analyte in the organic layer ( $C_o$ ) divided by the concentration of the analyte in the water layer ( $C_w$ ).

The amount of water should be the same before and after the extraction, because water is insoluble in 1-octanol. Therefore the two integral values we have obtained for the analyte are directly proportional to the concentration before and after the extraction.

Thus, the following equation can be used:

$$P_{ow} = \frac{C_o}{C_w} = \frac{I_w - I_{wo}}{I_w},$$

where  $I_w$  and  $I_{wo}$  are the peak values of the analyte in water and in water with 1-octanol added.

## Example $^1\text{H}$ NMR spectra

### Acetone




Figure 2. Spectrum of acetone in water (top) and with a layer of 1-octanol (bottom).

$$P_{ow} = C_o/C_w = \frac{0.197 - 0.124}{0.124} = 0.589$$

$$\log P_{ow} = -0.23 \text{ (literature value} = -0.24)^2$$

## Ethanol




Figure 3. Spectrum of ethanol in water (top) and with a layer of 1-octanol (bottom).

$$P_{ow} = C_o/C_w = \frac{0.183 - 0.121}{0.121} = 0.512$$

$$\log P_{ow} = -0.29 \text{ (literature value} = -0.30)^2$$

## THF




Figure 4. Spectrum of THF in water (top) and with a layer of 1-octanol (bottom).

$$P_{ow} = C_o/C_w = \frac{0.152 - 0.041}{0.041} = 2.707$$

$$\log P_{ow} = 0.43 \text{ (literature value} = 0.46)^2$$

## Tert-butanol




Figure 5. Spectrum of *tert*-butanol in water (top) and with a layer of 1-octanol (bottom).

$$P_{ow} = C_o/C_w = \frac{0.197 - 0.060}{0.060} = 2.28$$

$$\log P_{ow} = 0.36 \text{ (literature value} = 0.35)^2$$

## 2-Propanol




Figure 6. Spectrum of 2-propanol in water (top) and with a layer of 1-octanol (bottom).

$$P_{ow} = C_o/C_w = \frac{0.249 - 0.109}{0.109} = 1.28$$

$$\log P_{ow} = 0.11 \text{ (literature value} = 0.05)^2$$

## Acetonitrile




Figure 7. Spectrum of acetonitrile in water (top) and with a layer of 1-octanol (bottom).

$$P_{ow} = C_o/C_w = \frac{0.098 - 0.075}{0.075} = 0.307$$

$$\text{Log } P_{ow} = -0.51 \text{ (literature value} = -0.34)^2$$



# Spinsolve®



## References

- 1) Lipinski, C. A.; Lombardo, F.; Dominy, B. W.; Feeney, P. J., Experimental and computational approaches to estimate solubility and permeability in drug discovery and development settings. *Adv. Drug Deliv. Rev.* **1997**, 23 (1-3), 3-25.
- 2) Sangster, J., Octanol-water coefficients of simple organic compounds. *J. Phys. Chem. Ref. Data.* **1989**, 18, 1111-1227.

## CONTACT INFORMATION

For further information, please contact: [sales@magritek.com](mailto:sales@magritek.com)

### UNITED STATES

6440 Lusk Blvd (D108)  
San Diego, CA 92121, USA  
Tel: (855) 667-6835  
(855) NMR-MTEK

### GERMANY

Philipsstraße 8  
52068 Aachen, Germany  
Tel: +49 (241) 70525-6000  
Fax: +49 (241) 9278-6939

### NEW ZEALAND

32 Salamanca Road  
Wellington 6012, NZ  
Tel: +64 (4) 920-7671  
Fax: +64 (4) 471-4665

Or visit our website [www.magritek.com](http://www.magritek.com)